

## Peroxisome proliferator-activated receptor delta

Catalog No: #AP86491



Package Size: #AP86491-1 50ug #AP86491-2 100ug #AP86491-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Peroxisome proliferator-activated receptor delta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PPAR-delta,Nuclear hormone receptor 1,NUC1,Nuclear receptor subfamily 1 group C member 2,Peroxisome proliferator-activated receptor beta,PPAR-beta,Nr1c2,Pparb
Accession No.	Uniprot ID:P35396
Uniprot	P35396
GeneID	19015;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ppard

Note: This product is for in vitro research use only