

## Protein O-mannosyl-transferase 1

Catalog No: #AP86498



Package Size: #AP86498-1 50ug #AP86498-2 100ug #AP86498-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Protein O-mannosyl-transferase 1                          |
| Brief Description  | Recombinant Protein                                       |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE                              |
| Species Reactivity | Ms  |
| Other Names        | Dolichyl-phosphate-mannose--protein mannosyltransferase 1 |
| Accession No.      | Uniprot ID:Q8R2R1   |
| Uniprot            | Q8R2R1  |
| GeneID             | 99011;  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)    |

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Pomt1

Note: This product is for in vitro research use only