

Serine/threonine-protein phosphatase PP1-alpha catalytic subunit

Catalog No: #AP86500

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86500-1 50ug #AP86500-2 100ug #AP86500-3 1mg

Description

Product Name	Serine/threonine-protein phosphatase PP1-alpha catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	PP-1A
Accession No.	Uniprot ID:Q8WMS6
Uniprot	Q8WMS6
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PPP1CA

Note: This product is for in vitro research use only