

## Phospholipase D3

Catalog No: #AP86502



Package Size: #AP86502-1 50ug #AP86502-2 100ug #AP86502-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Phospholipase D3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PLD 3,Choline phosphatase 3,Phosphatidylcholine-hydrolyzing phospholipase D3,Schwannoma-associated protein 9,SAM-9,Sam9
Accession No.	Uniprot ID:O35405
Uniprot	O35405
GeneID	18807;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pld3

Note: This product is for in vitro research use only