

Myosin-2

Catalog No: #AP86503



Package Size: #AP86503-1 50ug #AP86503-2 100ug #AP86503-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Myosin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Myosin heavy chain 2,Myosin heavy chain 2a,MyHC-2a,Myosin heavy chain, skeletal muscle, adult 2
Accession No.	Uniprot ID:Q9TV63
Uniprot	Q9TV63
GeneID	397256;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:MYH2

Note: This product is for in vitro research use only