

## Platelet-derived growth factor receptor beta

Catalog No: #AP86506



Package Size: #AP86506-1 50ug #AP86506-2 100ug #AP86506-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Platelet-derived growth factor receptor beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PDGF-R-beta,Beta platelet-derived growth factor receptor,Beta-type platelet-derived growth factor receptor,CD140 antigen-like family member B,Platelet-derived growth factor receptor 1,PDGFR-1,CD140b,PDGFR,PDGFR1
Accession No.	Uniprot ID:P09619
Uniprot	P09619
GeneID	5159;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PDGFRB

Note: This product is for in vitro research use only