

# ProSAAS

Catalog No: #AP86507



Package Size: #AP86507-1 50ug #AP86507-2 100ug #AP86507-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	ProSAAS
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Proprotein convertase subtilisin/kexin type 1 inhibitor,Proprotein convertase 1 inhibitor,pro-SAAS
Accession No.	Uniprot ID:Q9QXU9
Uniprot	Q9QXU9
GeneID	246333;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pcsk1n

Note: This product is for in vitro research use only