

Protein disulfide-isomerase A3

Catalog No: #AP86509



Package Size: #AP86509-1 50ug #AP86509-2 100ug #AP86509-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein disulfide-isomerase A3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	58 kDa glucose-regulated protein,58 kDa microsomal protein,p58,Disulfide isomerase ER-60,Endoplasmic reticulum resident protein 57,ER protein 57,Endoplasmic reticulum resident protein 60,ER protein 60,HIP-70,Q-2,Erp60,Grp58
Accession No.	Uniprot ID:P11598
Uniprot	P11598
GeneID	29468;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdia3

Note: This product is for in vitro research use only