

## Proliferation-associated protein 2G4

Catalog No: #AP86511



Package Size: #AP86511-1 50ug #AP86511-2 100ug #AP86511-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Proliferation-associated protein 2G4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cell cycle protein p38-2G4 homolog,hG4-1,ErbB3-binding protein 1,EBP1
Accession No.	Uniprot ID:Q9UQ80
Uniprot	Q9UQ80
GeneID	5036;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PA2G4

Note: This product is for in vitro research use only