

Protein-arginine deiminase type-2

Catalog No: #AP86513



Package Size: #AP86513-1 50ug #AP86513-2 100ug #AP86513-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein-arginine deiminase type-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Peptidylarginine deiminase II,Protein-arginine deiminase type II,Pad2,Pdi,Pdi2
Accession No.	Uniprot ID:Q08642
Uniprot	Q08642
GeneID	18600;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Padi2

Note: This product is for in vitro research use only