

Prostate androgen-regulated mucin-like protein 1 homolog

Catalog No: #AP86515

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86515-1 50ug #AP86515-2 100ug #AP86515-3 1mg

Description

Product Name	Prostate androgen-regulated mucin-like protein 1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PARM-1,Castration-induced prostatic apoptosis-related protein 1,CIPAR-1,Prostatic androgen-repressed message 1 protein,Cipar1
Accession No.	Uniprot ID:Q6P9X9
Uniprot	Q6P9X9
GeneID	286894;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Parm1

Note: This product is for in vitro research use only