

Nucleolar and spindle-associated protein 1

Catalog No: #AP86517



Package Size: #AP86517-1 50ug #AP86517-2 100ug #AP86517-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nucleolar and spindle-associated protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	NuSAP,RCJMB04_15k16
Accession No.	Uniprot ID:Q5ZJU5
Uniprot	Q5ZJU5
GeneID	423213;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NUSAP1

Note: This product is for in vitro research use only