

Nuclear factor NF-kappa-B p100 subunit

Catalog No: #AP86522



Package Size: #AP86522-1 50ug #AP86522-2 100ug #AP86522-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Nuclear factor NF-kappa-B p100 subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	DNA-binding factor KBF2,Lyt-10,Nuclear factor NF-kappa-B p97 subunit,Nuclear factor of kappa light polypeptide gene enhancer in B-cells 2
Accession No.	Uniprot ID:P98150
Uniprot	P98150
GeneID	386574;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:NFKB2

Note: This product is for in vitro research use only