

Platelet-activating factor acetylhydrolase

Catalog No: #AP86526



Package Size: #AP86526-1 50ug #AP86526-2 100ug #AP86526-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-activating factor acetylhydrolase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PAF acetylhydrolase,1-alkyl-2-acetyl-glycerophosphocholine esterase,2-acetyl-1-alkyl-glycerophosphocholine esterase,LDL-associated phospholipase A2,LDL-PLA(2),PAF 2-acylhydrolase,Pafah
Accession No.	Uniprot ID:Q60963
Uniprot	Q60963
GeneID	27226;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pla2g7

Note: This product is for in vitro research use only