

Lysophosphatidylcholine acyltransferase 2

Catalog No: #AP86528

Package Size: #AP86528-1 50ug #AP86528-2 100ug #AP86528-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysophosphatidylcholine acyltransferase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LPC acyltransferase 2,1-acylglycerol-3-phosphate O-acyltransferase 11,1-AGP acyltransferase 11,1-acylglycerophosphocholine O-acyltransferase,1-alkylglycerophosphocholine O-acetyltransferase,Acetyl-CoA:lyso-platelet-activating factor acetyltransferase,Acetyl-CoA:lyso-PAF acetyltransferase,Acyltransferase-like 1,Lysophosphatidic acid acyltransferase alpha,LPAAT-alpha,AGPAT11,AYTL1
Accession No.	Uniprot ID:Q7L5N7
Uniprot	Q7L5N7
GeneID	54947;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LPCAT2

Note: This product is for in vitro research use only