

Phosphoenolpyruvate carboxykinase, cytosolic [GTP]

Catalog No: #AP86530



Package Size: #AP86530-1 50ug #AP86530-2 100ug #AP86530-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphoenolpyruvate carboxykinase, cytosolic [GTP]
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PEPCK-C,Pepck
Accession No.	Uniprot ID:Q9Z2V4
Uniprot	Q9Z2V4
GeneID	18534;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pck1

Note: This product is for in vitro research use only