

Retrotransposon-derived protein PEG10

Catalog No: #AP86532



Package Size: #AP86532-1 50ug #AP86532-2 100ug #AP86532-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Retrotransposon-derived protein PEG10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Embryonal carcinoma differentiation regulated protein, Mammalian retrotransposon-derived protein 2, Myelin expression factor 3, MyEF-3, Myelin expression factor 3-like protein 1, MEF3-like protein 1, Paternally expressed gene 10 protein, Retrotransposon gag domain-containing protein 3, Retrotransposon-derived gag-like polyprotein, Ty3/Gypsy-like protein, Edr, Mar2, Mart2, Mef311, Rgag3
Accession No.	Uniprot ID:Q7TN75
Uniprot	Q7TN75
GeneID	170676;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Peg10

Note: This product is for in vitro research use only