

Prostaglandin F2-alpha receptor

Catalog No: #AP86533



Package Size: #AP86533-1 50ug #AP86533-2 100ug #AP86533-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prostaglandin F2-alpha receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PGF receptor, Prostanoid FP receptor
Accession No.	Uniprot ID:P43117
Uniprot	P43117
GeneID	19220;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Ptgfr

Note: This product is for in vitro research use only