

PHD finger protein 1

Catalog No: #AP86536



Package Size: #AP86536-1 50ug #AP86536-2 100ug #AP86536-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PHD finger protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Protein PHF1, Polycomb-like protein 1, hPCI1, PCL1
Accession No.	Uniprot ID:O43189
Uniprot	O43189
GeneID	5252;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PHF1

Note: This product is for in vitro research use only