

POU domain, class 2, transcription factor 2

Catalog No: #AP86537



Package Size: #AP86537-1 50ug #AP86537-2 100ug #AP86537-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	POU domain, class 2, transcription factor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lymphoid-restricted immunoglobulin octamer-binding protein NF-A2, Octamer-binding protein 2, Oct-2, Octamer-binding transcription factor 2, OTF-2, OCT2, OTF2
Accession No.	Uniprot ID:P09086
Uniprot	P09086
GeneID	5452;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: POU2F2

Note: This product is for in vitro research use only