

Platelet-derived growth factor receptor alpha

Catalog No: #AP86538



Package Size: #AP86538-1 50ug #AP86538-2 100ug #AP86538-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor receptor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	PDGF-R-alpha, Alpha platelet-derived growth factor receptor, Alpha-type platelet-derived growth factor receptor
Accession No.	Uniprot ID:Q9PUF6
Uniprot	Q9PUF6
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PDGFRA

Note: This product is for in vitro research use only