

Platelet-derived growth factor receptor alpha

Catalog No: #AP86539



Package Size: #AP86539-1 50ug #AP86539-2 100ug #AP86539-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor receptor alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PDGF-R-alpha,Alpha platelet-derived growth factor receptor,Alpha-type platelet-derived growth factor receptor,CD140 antigen-like family member A,Platelet-derived growth factor alpha receptor,CD140a
Accession No.	Uniprot ID:P26618
Uniprot	P26618
GeneID	18595;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdgfra

Note: This product is for in vitro research use only