

Peroxidasin homolog

Catalog No: #AP86543

Package Size: #AP86543-1 50ug #AP86543-2 100ug #AP86543-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxidasin homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Melanoma-associated antigen MG50,Vascular peroxidase 1,p53-responsive gene 2 protein,KIAA0230,MG50,PRG2,VPO,VPO1
Accession No.	Uniprot ID:Q92626
Uniprot	Q92626
GeneID	7837;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PXDN

Note: This product is for in vitro research use only