

# Platelet factor 4

Catalog No: #AP86544



Package Size: #AP86544-1 50ug #AP86544-2 100ug #AP86544-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Platelet factor 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PF-4,C-X-C motif chemokine 4,Iroplact,Oncostatin-A,CXCL4,SCYB4
Accession No.	Uniprot ID:P02776
Uniprot	P02776
GeneID	5196;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PF4

Note: This product is for in vitro research use only