

# Serum paraoxonase/arylesterase 2

Catalog No: #AP86546



Package Size: #AP86546-1 50ug #AP86546-2 100ug #AP86546-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Serum paraoxonase/arylesterase 2                                      |
| Brief Description  | Recombinant Protein   |
| Host Species       | E.coli  |
| Purification       | Greater than 90% by SDS-PAGE  |
| Species Reactivity | Dog   |
| Other Names        | PON 2,Aromatic esterase 2,A-esterase 2,Serum aryldialkylphosphatase 2 |
| Accession No.      | Uniprot ID:P54832   |
| Uniprot            | P54832  |
| Tag Info           | His   |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)                |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PON2

Note: This product is for in vitro research use only