

Oncoprotein-induced transcript 3 protein

Catalog No: #AP86547



Package Size: #AP86547-1 50ug #AP86547-2 100ug #AP86547-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oncoprotein-induced transcript 3 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Liver-specific zona pellucida domain-containing protein,LZP,UNQ826/PRO1753
Accession No.	Uniprot ID:Q8WWZ8
Uniprot	Q8WWZ8
GeneID	170392;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:OIT3

Note: This product is for in vitro research use only