

Protein kinase C-binding protein NELL2

Catalog No: #AP86548



Package Size: #AP86548-1 50ug #AP86548-2 100ug #AP86548-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein kinase C-binding protein NELL2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	MEL91 protein, NEL-like protein 2, Mel91
Accession No.	Uniprot ID:Q61220
Uniprot	Q61220
GeneID	54003;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Nell2

Note: This product is for in vitro research use only