

Proteinase-activated receptor 2

Catalog No: #AP86550



Package Size: #AP86550-1 50ug #AP86550-2 100ug #AP86550-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proteinase-activated receptor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PAR-2,Coagulation factor II receptor-like 1,G-protein coupled receptor 11,Thrombin receptor-like 1,Gpcr11,Gpr11,Par2
Accession No.	Uniprot ID:P55086
Uniprot	P55086
GeneID	14063;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:F2r11

Note: This product is for in vitro research use only