

# High affinity cGMP-specific 3',5'-cyclic phosphodiesterase 9A



Catalog No: #AP86552

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP86552-1 50ug #AP86552-2 100ug #AP86552-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	High affinity cGMP-specific 3',5'-cyclic phosphodiesterase 9A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pde8b
Accession No.	Uniprot ID:O70628
Uniprot	O70628
GeneID	18585;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Pde9a

Note: This product is for in vitro research use only