

Platelet-derived growth factor C

Catalog No: #AP86554



Package Size: #AP86554-1 50ug #AP86554-2 100ug #AP86554-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PDGF-C, Fallotein, Spinal cord-derived growth factor, SCDGF, VEGF-E, SCDGF, UNQ174/PRO200
Accession No.	Uniprot ID:Q9NRA1
Uniprot	Q9NRA1
GeneID	56034;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: PDGFC

Note: This product is for in vitro research use only