

Prolyl 4-hydroxylase subunit alpha-2

Catalog No: #AP86564



Package Size: #AP86564-1 50ug #AP86564-2 100ug #AP86564-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prolyl 4-hydroxylase subunit alpha-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	4-PH alpha-2,Procollagen-proline,2-oxoglutarate-4-dioxygenase subunit alpha-2
Accession No.	Uniprot ID:Q60716
Uniprot	Q60716
GeneID	18452;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: P4ha2

Note: This product is for in vitro research use only