

Cellular tumor antigen p53

Catalog No: #AP86565



Package Size: #AP86565-1 50ug #AP86565-2 100ug #AP86565-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cellular tumor antigen p53
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Tumor suppressor p53,P53
Accession No.	Uniprot ID:Q9TUB2
Uniprot	Q9TUB2
GeneID	397276;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:TP53

Note: This product is for in vitro research use only