

Serine/threonine-protein kinase PLK2

Catalog No: #AP86566



Package Size: #AP86566-1 50ug #AP86566-2 100ug #AP86566-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase PLK2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Polo-like kinase 2,PLK-2,Serine/threonine-protein kinase SNK,Serum-inducible kinase,Snk
Accession No.	Uniprot ID:P53351
Uniprot	P53351
GeneID	20620;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Plk2

Note: This product is for in vitro research use only