

Platelet-derived growth factor subunit B

Catalog No: #AP86568



Package Size: #AP86568-1 50ug #AP86568-2 100ug #AP86568-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor subunit B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	PDGF subunit B,PDGF-2,Platelet-derived growth factor B chain,Platelet-derived growth factor beta polypeptide
Accession No.	Uniprot ID:Q05028
Uniprot	Q05028
GeneID	24628;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pdgfb

Note: This product is for in vitro research use only