

Platelet-derived growth factor D

Catalog No: #AP86569



Package Size: #AP86569-1 50ug #AP86569-2 100ug #AP86569-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor D
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PDGF-D, Spinal cord-derived growth factor B, SCDGF-B, Scdgb
Accession No.	Uniprot ID:Q92517
Uniprot	Q92517
GeneID	71785;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pdgfd

Note: This product is for in vitro research use only