

Platelet-derived growth factor subunit A

Catalog No: #AP86574



Package Size: #AP86574-1 50ug #AP86574-2 100ug #AP86574-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Platelet-derived growth factor subunit A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PDGF subunit A,Platelet-derived growth factor A chain,Platelet-derived growth factor alpha polypeptide
Accession No.	Uniprot ID:Q2KJ15
Uniprot	Q2KJ15
GeneID	505908;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PDGFA

Note: This product is for in vitro research use only