

CDP-diacylglycerol--glycerol-3-phosphate 3-phosphatidyltransferase, mitochondrial

Catalog No: #AP86577

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP86577-1 50ug #AP86577-2 100ug #AP86577-3 1mg

Description

Product Name	CDP-diacylglycerol--glycerol-3-phosphate 3-phosphatidyltransferase, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Phosphatidylglycerophosphate synthase 1,PGP synthase 1,Silencer-associated factor,Saf
Accession No.	Uniprot ID:Q8BHF7
Uniprot	Q8BHF7
GeneID	74451;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pgs1

Note: This product is for in vitro research use only