

Phosphoenolpyruvate carboxykinase [GTP], mitochondrial

Catalog No: #AP86580

Package Size: #AP86580-1 50ug #AP86580-2 100ug #AP86580-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Phosphoenolpyruvate carboxykinase [GTP], mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	PEPCK-M
Accession No.	Uniprot ID:P21642
Uniprot	P21642
GeneID	396457;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PCK2

Note: This product is for in vitro research use only