

Peroxiredoxin-5, mitochondrial

Catalog No: #AP86583



Package Size: #AP86583-1 50ug #AP86583-2 100ug #AP86583-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Peroxiredoxin-5, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Antioxidant enzyme B166,AOEB166,PLP,Peroxiredoxin V,Prx-V,Peroxisomal antioxidant enzyme,Thioredoxin peroxidase PMP20,Thioredoxin reductase
Accession No.	Uniprot ID:Q9R063
Uniprot	Q9R063
GeneID	113898;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prdx5

Note: This product is for in vitro research use only