

Oxidative stress-induced growth inhibitor 2

Catalog No: #AP86585



Package Size: #AP86585-1 50ug #AP86585-2 100ug #AP86585-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Oxidative stress-induced growth inhibitor 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	hT41,C8orf1
Accession No.	Uniprot ID:Q9Y236
Uniprot	Q9Y236
GeneID	734;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: OSGIN2

Note: This product is for in vitro research use only