

Natural resistance-associated macrophage protein 2

Catalog No: #AP86587



Package Size: #AP86587-1 50ug #AP86587-2 100ug #AP86587-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Natural resistance-associated macrophage protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	NRAMP 2,Divalent cation transporter 1,Divalent metal transporter 1,DMT-1,Solute carrier family 11 member 2,Dct1,Dmt1,Nramp2
Accession No.	Uniprot ID:P49282
Uniprot	P49282
GeneID	18174;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Slc11a2

Note: This product is for in vitro research use only