

Lysosomal protective protein

Catalog No: #AP86588

Package Size: #AP86588-1 50ug #AP86588-2 100ug #AP86588-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysosomal protective protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Cathepsin A,PPGB
Accession No.	Uniprot ID:Q3MI05
Uniprot	Q3MI05
GeneID	518169;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CTSA

Note: This product is for in vitro research use only