

# Peroxisome proliferator-activated receptor gamma coactivator-related protein 1

Catalog No: #AP86589

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP86589-1 50ug #AP86589-2 100ug #AP86589-3 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Peroxisome proliferator-activated receptor gamma coactivator-related protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PGC-1-related coactivator,PRC,KIAA0595
Accession No.	Uniprot ID:Q5VV67
Uniprot	Q5VV67
GeneID	23082;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PPRC1

Note: This product is for in vitro research use only