

Serine/threonine-protein phosphatase PP1-alpha catalytic subunit

Catalog No: #AP86590

Package Size: #AP86590-1 50ug #AP86590-2 100ug #AP86590-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein phosphatase PP1-alpha catalytic subunit
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PP-1A,Ppp1a
Accession No.	Uniprot ID:P62137
Uniprot	P62137
GeneID	19045;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppp1ca

Note: This product is for in vitro research use only