

Placental protein 13-like

Catalog No: #AP86591



Package Size: #AP86591-1 50ug #AP86591-2 100ug #AP86591-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Placental protein 13-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Charcot-Leyden crystal protein 2,CLC2,Galectin-14,Gal-14,PPL13
Accession No.	Uniprot ID:Q8TCE9
Uniprot	Q8TCE9
GeneID	56891;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LGALS14

Note: This product is for in vitro research use only