

Oxytocin receptor

Catalog No: #AP86593



Package Size: #AP86593-1 50ug #AP86593-2 100ug #AP86593-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Oxytocin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	OT-R,Otr
Accession No.	Uniprot ID:P70536
Uniprot	P70536
GeneID	25342;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Oxtr

Note: This product is for in vitro research use only