

Prolyl 3-hydroxylase 1

Catalog No: #AP86594



Package Size: #AP86594-1 50ug #AP86594-2 100ug #AP86594-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Prolyl 3-hydroxylase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Growth suppressor 1,Leucine- and proline-enriched proteoglycan 1,Leprecan-1,Gros1,Lepre1
Accession No.	Uniprot ID:Q3V1T4
Uniprot	Q3V1T4
GeneID	56401;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:P3h1

Note: This product is for in vitro research use only