

# Pituitary adenylate cyclase-activating polypeptide type I receptor

Catalog No: #AP86595

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP86595-1 50ug #AP86595-2 100ug #AP86595-3 1mg

## Description

Product Name	Pituitary adenylate cyclase-activating polypeptide type I receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	PACAP type I receptor
Accession No.	Uniprot ID:Q29627
Uniprot	Q29627
GeneID	319095;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ADCYAP1R1

Note: This product is for in vitro research use only