

PCNA-associated factor

Catalog No: #AP86596



Package Size: #AP86596-1 50ug #AP86596-2 100ug #AP86596-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PCNA-associated factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Hepatitis C virus NS5A-transactivated protein 9,HCV NS5A-transactivated protein 9,Overexpressed in anaplastic thyroid carcinoma 1,OEATC-1,p15PAF,KIAA0101,NS5ATP9,L5
Accession No.	Uniprot ID:Q15004
Uniprot	Q15004
GeneID	9768;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PCLAF

Note: This product is for in vitro research use only