

Calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1C

Catalog No: #AP86599

Orders: order@signalwayantibody.com

Package Size: #AP86599-1 50ug #AP86599-2 100ug #AP86599-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Calcium/calmodulin-dependent 3',5'-cyclic nucleotide phosphodiesterase 1C
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cam-PDE 1C
Accession No.	Uniprot ID:Q64338
Uniprot	Q64338
GeneID	18575;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Pde1c

Note: This product is for in vitro research use only